

## How it works

Simply back your wheelchair up to the **SPEEDYLIFT**, push the 'IN' button and the rest is done for you. The **SPEEDYLIFT** picks up your chair and tucks it into place behind the driver's seat.

## Features

The compact design takes up very little space. In most cases, only one seat is removed from the inside of your vehicle. That means you can keep the use of up to 6 seats!

Seeing the intriguing **SPEEDYLIFT** in action, you will understand where it gets its name from: A full cycle takes less than 30 seconds.

## Specifications

- Maximum capacity: 95 lbs
- An automatic locking mechanism fastens your wheelchair in place as it enters the vehicle and releases on the way out.


## Installation

The **SPEEDYLIFT** is easily installed into the second row seat mounts with no structural modifications to your vehicle. It is also easily removed to return the vehicle to its original condition without affecting resale value.

## Options

- Although it was designed for rigid frame manual wheelchairs, an optional adaptor kit is available for manual folding frame wheelchairs.
- Can be paired with the **XL-SEAT** or the **XL-BOARD**





# SPEEDY-LIFT

Fast wheelchair lift

The **SPEEDY-LIFT** is designed to fit adjacent to the sliding door in most minivans as well as the Mazda 5 on the driver or passenger side.

Because of its unparalleled speed, precision and ease, the **SPEEDY-LIFT** is an ideal solution for someone looking for simplified vehicle access.

## Safety

We take our customers' safety seriously. We are dedicated to manufacturing quality mobility equipment. Our products undergo a strict certification process. Certification of our products involves different tests such as: crash testing, cycle testing and weight capacity testing. We are proud to offer a line of products that exceed the industry's standards.



## Vehicle Applications

- Buick - Terraza / Silhouette
- Chevrolet - Uplander / Venture
- Chrysler - Town & Country
- Dodge - Caravan
- Ford - Windstar / Freestar
- Honda - Odyssey
- Hyundai - Entourage
- Kia - Sedona
- Mazda - Mazda 5
- Pontiac - Transport / Montana / Montana SV6
- Saturn - Relay
- Toyota - Sienna
- Volkswagen - Routan

*Due to the diversity of mobility solutions available to you, please contact your local mobility dealer for a personalized evaluation and overview of options tailored to your needs.*

Our products qualify as adaptive equipment and are eligible for reimbursement under the terms of Mobility programs offered by various vehicle manufacturers.

Dealer :

## ADAPT SOLUTIONS

QUALITY MOBILITY EQUIPMENT

145 DAMASE-BRETON | SAINT-LAMBERT, QC, G0S 2W0

P 866 641 0419 | F 418 889 9838 | O 418 889 0419

[www.adaptsolutions.ca](http://www.adaptsolutions.ca)



creation : [www.noirdstudio.com](http://www.noirdstudio.com)

[www.adaptsolutions.ca](http://www.adaptsolutions.ca)